


 COROENT


## CoRoent Small System



Small  
14 x 11mm



Small Lordotic Plus  
15 x 13mm



Small Lordotic  
14 x 11mm



The Total Solution  
for Cervical  
Interbody Fusion

## CoRoent Small Interbody Portfolio



### Multiple Choices, One Set

#### CoRoent Small, Small Lordotic, and Small Lordotic Plus

The complete line of CoRoent Small Cervical interbodies from NuVasive is precision-machined from PEEK, which balances structural integrity with elasticity properties closely matched to those of cortical bone. Radiolucency, along with titanium visualization markers, provide an ideal combination for proper alignment and fusion assessment. All of the CoRoent implant features are designed with patient safety and improved outcomes in mind.



#### CoRoent Small

- 14 x 11mm footprint
- Parallel
- Available in 5-12mm heights





#### CoRoent Small Lordotic

- 14 x 11mm footprint
- 5° lordotic taper
- Available in 5-12mm heights



#### CoRoent Small Lordotic Plus


- 15 x 13mm footprint
- 5° lordotic taper
- Available in 5-12mm heights



### CoRoent SLP Features

The new CoRoent SLP from NuVasive features an extended 15 x 13mm footprint, creating additional stability while maintaining a generous graft window for an optimal fusion environment.\*

- 1 Extended 15 x 13mm footprint for a new level of stability and subsidence resistance.
- 2 Market-leading central aperture size provides the ideal environment for fast and complete fusion.
- 3 Lordotic taper and directional biting teeth match disc space angulation and prevent migration.



\*CoRoent SLP's 15 x 13mm footprint was specifically designed to resist subsidence by resting on the denser bone of the posterior and lateral areas of the vertebral body.

# CoRoent Small System Instrumentation and Tray

The CoRoent Small System was designed to give the surgeon intraoperative choices to match any patient anatomy or surgical situation. The CoRoent Small System includes CoRoent Small, Small Lordotic, and Small Lordotic Plus implants, as well as matching trials and full instrumentation.



## Single Level Implant Access


CoRoent Small, Small Lordotic, and Small Lordotic Plus implants are all accessible on the first level of the CoRoent Small System, giving the surgeon maximal access to a wide variety of options to accommodate any surgical situation.

## Color-Coded, Low-Profile Trials

A full selection of trials is included in the CoRoent Small System to complement all implant offerings. Trials are color-coded and low-profile for maximum visibility and ease of use.



## CoRoent Small System



## Additional Instrumentation

Also included in the CoRoent Small System are: implant tamp, shortened inserter, standard inserter, mallet, and core material insertion block. Maximizing interbody and instrument choices allows the surgeon to select the proper tool for varying surgical situations.

A cervical rasp tray is also available with both parallel and lordotic rasps to match all the NuVasive cervical offerings.

# CoRoent Small System Implants

## COROENT SMALL



| DESCRIPTION                  | CATALOG # |
|------------------------------|-----------|
| 5 x 14 x 11mm CoRoent Small  | 6805115   |
| 6 x 14 x 11mm CoRoent Small  | 6806115   |
| 7 x 14 x 11mm CoRoent Small  | 6807115   |
| 8 x 14 x 11mm CoRoent Small  | 6808115   |
| 9 x 14 x 11mm CoRoent Small  | 6809115   |
| 10 x 14 x 11mm CoRoent Small | 6810115   |
| 11 x 14 x 11mm CoRoent Small | 6811115   |
| 12 x 14 x 11mm CoRoent Small | 6812115   |

## COROENT SMALL LORDOTIC



| DESCRIPTION                           | CATALOG # |
|---------------------------------------|-----------|
| 5 x 14 x 11mm CoRoent Small Lordotic  | 6805714   |
| 6 x 14 x 11mm CoRoent Small Lordotic  | 6806714   |
| 7 x 14 x 11mm CoRoent Small Lordotic  | 6807714   |
| 8 x 14 x 11mm CoRoent Small Lordotic  | 6808714   |
| 9 x 14 x 11mm CoRoent Small Lordotic  | 6809714   |
| 10 x 14 x 11mm CoRoent Small Lordotic | 6810714   |
| 11 x 14 x 11mm CoRoent Small Lordotic | 6811714   |
| 12 x 14 x 11mm CoRoent Small Lordotic | 6812714   |

## COROENT SMALL LORDOTIC PLUS



| DESCRIPTION                | CATALOG # |
|----------------------------|-----------|
| 5 x 15 x 13mm CoRoent SLP  | 6805815   |
| 6 x 15 x 13mm CoRoent SLP  | 6806815   |
| 7 x 15 x 13mm CoRoent SLP  | 6807815   |
| 8 x 15 x 13mm CoRoent SLP  | 6808815   |
| 9 x 15 x 13mm CoRoent SLP  | 6809815   |
| 10 x 15 x 13mm CoRoent SLP | 6810815   |
| 11 x 15 x 13mm CoRoent SLP | 6811815   |
| 12 x 15 x 13mm CoRoent SLP | 6812815   |

**PLEASE NOTE: The NuVasive CoRoent Small family of interbodies is now cleared for use in the cervical spine. Take confidence in the fact that our CoRoent Small implants were specifically designed to complement the anatomical requirements of the cervical spine.**



To order, please contact your NuVasive Sales Consultant or Customer Service Representative today at:

**NuVasive, Inc.** 7475 Lusk Blvd., San Diego, CA 92121 USA • phone: 800-475-9131 fax: 800-475-9134



**NuVasive Netherlands B.V.** Jachthavenweg 109A, 1081 KM Amsterdam, The Netherlands • phone: +31 20 72 33 000

[nuvasive.com](http://nuvasive.com)