



CoRoent


CoRoent XL & XLR



XL



XL-Thoracic



XLR

Sizable
Advantages.


Sized for Stability – Shaped for Precise Fit – Designed for Simplicity.

CoRoent Features

- Large, anatomical shape for maximum surface area and structural stability
- Large apertures for bony through-growth
- Radiolucent PEEK-OPTIMA material provides optimal stiffness compatibility with surrounding bone
- Titanium bars enable easy-to-interpret placement and orientation verification



XL




- Designed for the unique eXtreme Lateral approach
- Multiple-length options to ensure optimal apophyseal support, reducing the chance of subsidence
- Available in standard and lordotic profiles to induce proper sagittal alignment




XL-Thoracic



- Smaller sizes to accommodate thoracic anatomy using the eXtreme Lateral approach
- Specially designed instruments for eXtreme Lateral thoracic surgery
- Single aperture for maximum through-growth
- Rounded anterior corners for ideal anatomical disc space coverage



XLR



- Designed for anterior or anterolateral approach
- Dual inserter attachment points for approach-specific placement
- Two footprint sizes for optimal stability and surface contact
- 8° and 15° lordosis options
- Specialized instruments and a unique delivery system to facilitate placement

CoRoent XL-Standard

XL-STANDARD 40, 45mm	CAT #	XL-STANDARD 50, 55mm	CAT #
8 x 18 x 40mm	6908840	8 x 18 x 50mm	6908850
10 x 18 x 40mm	6910840	10 x 18 x 50mm	6910850
12 x 18 x 40mm	6912840	12 x 18 x 50mm	6912850
14 x 18 x 40mm	6914840	14 x 18 x 50mm	6914850
16 x 18 x 40mm	6916840	16 x 18 x 50mm	6916850
8 x 18 x 45mm	6908845	8 x 18 x 55mm	6908855
10 x 18 x 45mm	6910845	10 x 18 x 55mm	6910855
12 x 18 x 45mm	6912845	12 x 18 x 55mm	6912855
14 x 18 x 45mm	6914845	14 x 18 x 55mm	6914855
16 x 18 x 45mm	6916845	16 x 18 x 55mm	6916855

CoRoent XL-Lordotic

XL-LORDOTIC 40, 45mm	CAT #	XL-LORDOTIC 50, 55mm	CAT #
8 x 18 x 40mm 10°	6980840	8 x 18 x 50mm 10°	6980850
10 x 18 x 40mm 10°	6981040	10 x 18 x 50mm 10°	6981050
12 x 18 x 40mm 10°	6981240	12 x 18 x 50mm 10°	6981250
14 x 18 x 40mm 10°	6981440	14 x 18 x 50mm 10°	6981450
16 x 18 x 40mm 10°	6981640	16 x 18 x 50mm 10°	6981650
8 x 18 x 45mm 10°	6980845	8 x 18 x 55mm 10°	6980855
10 x 18 x 45mm 10°	6981045	10 x 18 x 55mm 10°	6981055
12 x 18 x 45mm 10°	6981245	12 x 18 x 55mm 10°	6981255
14 x 18 x 45mm 10°	6981445	14 x 18 x 55mm 10°	6981455
16 x 18 x 45mm 10°	6981645	16 x 18 x 55mm 10°	6981655

CoRoent XL-Thoracic

XL-THORACIC 25, 30, 35, 40mm	CAT #
6 x 16 x 25mm	6990625
8 x 16 x 25mm	6990825
10 x 16 x 25mm	6991025
12 x 16 x 25mm	6991225
6 x 16 x 30mm	6990630
8 x 16 x 30mm	6990830
10 x 16 x 30mm	6991030
12 x 16 x 30mm	6991230
6 x 16 x 35mm	6990635
8 x 16 x 35mm	6990835
10 x 16 x 35mm	6991035
12 x 16 x 35mm	6991235
6 x 16 x 40mm	6990640
8 x 16 x 40mm	6990840
10 x 16 x 40mm	6991040
12 x 16 x 40mm	6991240

CoRoent XLR

XLR – 34 x 24mm	CAT #
10 x 34 x 24mm 8°	6921044
12 x 34 x 24mm 8°	6921244
14 x 34 x 24mm 8°	6921444
16 x 34 x 24mm 8°	6921644
18 x 34 x 24mm 8°	6921844
20 x 34 x 24mm 8°	6922044
XLR – 38 x 28mm	CAT #
10 x 38 x 28mm 8°	6921088
12 x 38 x 28mm 8°	6921288
14 x 38 x 28mm 8°	6921488
16 x 38 x 28mm 8°	6921688
18 x 38 x 28mm 8°	6921888
20 x 38 x 28mm 8°	6922088
12 x 38 x 28mm 15°	6931288
14 x 38 x 28mm 15°	6931488
16 x 38 x 28mm 15°	6931688
18 x 38 x 28mm 15°	6931888
20 x 38 x 28mm 15°	6932088



To order, please contact your NuVasive Sales Consultant or Customer Service Representative today at:

NuVasive, Inc. 7475 Lusk Blvd., San Diego, CA 92121 USA • phone: 800-475-9131 fax: 800-475-9134

NuVasive Netherlands B.V. Jachthavenweg 109A, 1081 KM Amsterdam, The Netherlands • phone: +31 20 72 33 000

nuvasive.com

9004255 E

